

Determining the Concentration of a Solution: Beer's Law

The primary objective of this experiment is to determine the concentration of an unknown cobalt (II) chloride, CoCl_2 , solution. You will use a Vernier Spectrometer or a Vernier SpectroVis to measure the concentration of each solution. You will first measure the absorbance of a standard solution over the visible light spectrum and select the wavelength of maximum absorbance (λ_{max}). A higher concentration of the solution absorbs more light (and transmits less) than a solution of lower concentration.

You will prepare five cobalt (II) chloride solutions of known concentration (standard solutions) and measure the absorbance of the standard solutions at the λ_{max} . The graph of absorbance vs. concentration for the standard solutions will describe a direct relationship, known as *Beer's law*.

Finally, you will determine the concentration of an unknown CoCl_2 solution by measuring its absorbance with a spectrometer and using the best-fit line equation of the Beer's law curve to calculate the unknown's concentration.

OBJECTIVES

In this experiment, you will

- Measure and analyze the visible light absorbance spectrum of a standard cobalt (II) chloride solution to determine the maximum wavelength of absorbance.
- Prepare and test the absorbance of five standard cobalt (II) chloride solutions.
- Calculate a standard curve from the test results of the standard solutions.
- Test the absorbance of a cobalt (II) chloride solution of unknown molar concentration.


MATERIALS

| | |
|---|--|
| Vernier Spectrometer or SpectroVis | 0.20 M cobalt (II) chloride, CoCl_2 , solution |
| computer | cobalt (II) chloride, CoCl_2 , unknown solution |
| Vernier Logger <i>Pro</i> 3 software | distilled water |
| pipet pump or pipet bulb | test tube rack |
| one cuvette | stirring rod |
| two 10 mL pipets or graduated cylinders | tissues (preferably lint-free) |
| five 20 × 150 mm test tubes | two 100 mL beakers |


PROCEDURE

1. Obtain and wear goggles.
2. Obtain ~30 mL of 0.20 M CoCl₂ solution and distilled water.
3. Label four clean, dry, test tubes 1–4. Use pipets to prepare the standard solutions according to the chart below. Thoroughly mix each solution with a stirring rod. Clean and dry the stirring rod between uses.

| Test Tube number | 0.20 M CoCl ₂ (mL) | Distilled H ₂ O (mL) | Concentration (M) |
|------------------|-------------------------------|---------------------------------|-------------------|
| 1 | 10 | 0 | 0.20 |
| 2 | 8 | 2 | 0.16 |
| 3 | 6 | 4 | 0.12 |
| 4 | 4 | 6 | 0.08 |

4. Use a USB cable to connect a spectrometer to the computer.
5. Start the Logger Pro 3 (version 3.6 or newer) program on your computer.
6. Calibrate the spectrometer.
 - a. Prepare a *blank* by filling an empty cuvette $\frac{3}{4}$ full with distilled water. Place the blank cuvette in the spectrometer.
 - b. Select Calibrate ► Spectrometer from the Experiment menu. The calibration dialog box will display the message: “Waiting ...seconds for lamp to warm up.” The minimum warm up time is 90 seconds. Follow the instructions in the dialog box to complete the calibration. Click .
7. Determine the maximum wavelength for CoCl₂ (aq) and set up the data collection mode.
 - a. Empty the blank cuvette and rinse it twice with small amounts of the 0.20 M CoCl₂ solution in Test Tube 1. Fill the cuvette $\frac{3}{4}$ full with the 0.20 M CoCl₂ solution and place it in the spectrometer.
 - b. Click . A full spectrum graph of the CoCl₂ solution will be displayed. Note that one area of the graph contains a peak absorbance. Click to complete the analysis.
 - c. To save your graph of absorbance vs. wavelength, select Store Latest Run from the Experiment menu.
 - d. Click the Configure Spectrometer Data Collection icon, , on the toolbar. A dialog box will appear.
 - e. Select Abs vs. Concentration under Set Collection Mode. The wavelength of peak absorbance (λ max) will be automatically selected. If you wish to select a new wavelength, click on the graph or check the box next to the desired wavelength. Click to proceed.
8. Collect absorbance-concentration data for the five standard solutions.
 - a. Leave the cuvette in the spectrometer. Click . When the absorbance reading stabilizes, click . Enter “0.20” as the concentration of the solution and click .
 - b. Discard the cuvette contents as directed. Using the solution in Test Tube 2, rinse and fill the cuvette $\frac{3}{4}$ full. Wipe the cuvette and place it in the spectrometer. When the absorbance reading stabilizes, click . Enter “0.16” as the concentration.
 - c. Repeat Step 8b for the remaining test tubes of the standard CoCl₂ solution. When you have finished testing the standard solutions, click .

Determining the Concentration of a Solution: Beer's Law

9. To determine the best-fit line equation for the CoCl_2 standard solutions, click the linear fit button, , on the toolbar. Write down the equation for the standard solutions in your data table or lab book.
10. Determine the concentration of the unknown CoCl_2 solution.
 - a. Obtain about 5 mL of the *unknown* CoCl_2 solution.
 - b. Rinse the cuvette twice with the unknown solution and fill it about $\frac{3}{4}$ full. Wipe the outside of the cuvette and place it into the spectrometer.
 - c. Select Interpolation Calculator from the Analyze menu. A dialog box will appear that displays the concentration of your unknown at the measured absorbance.
 - d. Click . Write down the absorbance and the concentration of the unknown in your data table or lab book.
 - e. Dispose of any of the remaining solutions as directed.
11. (optional) Print a copy of your graphs and/or data table.
12. (optional) Select Save As... from the File menu and save your experiment file.

DATA TABLE

| Test Tube | Concentration (mol/L) | Absorbance |
|-----------|-----------------------|------------|
| 1 | 0.20 | |
| 2 | 0.16 | |
| 3 | 0.12 | |
| 4 | 0.08 | |
| 5 | Unknown number ____ | |

DATA ANALYSIS

Your lab report should contain the following information:

- a graph showing the data and linear-regression equation for the standard solutions
- the molar concentration of your unknown CoCl_2 solution
- an explanation of your calculation of the unknown's concentration